



## THE CHOICE IS YOURS



BEGO SEMADOS® S/SC/SCX IMPLANTS  
QUICK & FUNCTIONAL · UNIVERSAL & SUCCESSFUL  
RELIABLE & PROVEN · VARIED & SOPHISTICATED  
100% GERMAN QUALITY PRODUCT  
SIMPLE & INTUITIVE · QUICK & FUNCTIONAL  
UNIVERSAL & SUCCESSFUL · RELIABLE & PROVEN  
VARIED & SOPHISTICATED · SIMPLE & INTUITIVE

### Your tried and true system – familiar yet new

**Tailored to your patient down to the finest detail.** Various implant systems, supplemented by individual CAD/CAM fabricated prosthetic components and surgical guides as well as the BEGO Multi<sup>Plus</sup> concept – a system that provides patients with a perfectly fitting prosthesis in just one day. The exceptionally high technical standards ensure the highest level of satisfaction – for both dental implantologists and patients.

#### Your benefits

- S implants – increase in the primary stability thanks to increasing core diameter with decreasing thread depth
- SC/SCX implants – preservation of the crestal bone thanks to the integrated Platform Switch
- Efficient – just one surgical tray needed
- Well thought out – cover screw included
- Freedom of choice based on preference and indication – machined implant shoulder (S or SC implants) or micro-structured shoulder (SCX implants)
- Optional – guided/template-guided preparation of BEGO Semados® S/SC/SCX implants Ø 3.25–4.5 using the BEGO Guide System
- For the most stringent quality requirements – 100% German development and production

#### Technical details

- **Platform Match** (S implants) or **Platform Switch** (SC/SCX implants)
- Gap-free fit of the abutments thanks to the internal taper
- Self-tapping thread
- Manufactured from grade 4 pure titanium with ultrapure homogeneous TiPure<sup>Plus</sup> surface
- Diameters: 3.25 · 3.75 · 4.1 · 4.5 · 5.5 mm
- Lengths: 7 · 8.5 · 10 · 11.5 · 13 · 15 mm
- Diameter 5.5 mm in length 7 mm if there is not much vertical bone height available (SC/SCX only)

# PRODUCT OVERVIEW

Distributor contact data

BEGO Semados® S Implants			BEGO Semados® SC Implants			BEGO Semados® SCX Implants		
Art. no	Colour code	Product	Art. no	Colour code	Product	Art. no	Colour code	Product
<b>BEGO Semados® S Implants 3.25</b>			<b>BEGO Semados® SC Implants 3.25</b>			<b>BEGO Semados® SCX Implants 3.25</b>		
57500	■	S 3.25 L8.5	58160	■	SC 3.25 L8.5	58260	■	SCX 3.25 L8.5
57501	■	S 3.25 L10	58161	■	SC 3.25 L10	58261	■	SCX 3.25 L10
57502	■	S 3.25 L11.5	58162	■	SC 3.25 L11.5	58262	■	SCX 3.25 L11.5
57503	■	S 3.25 L13	58163	■	SC 3.25 L13	58263	■	SCX 3.25 L13
			58164	■	SC 3.25 L15	58264	■	SCX 3.25 L15
<b>BEGO Semados® S Implants 3.75</b>			<b>BEGO Semados® SC Implants 3.75</b>			<b>BEGO Semados® SCX Implants 3.75</b>		
57601	■	S 3.75 L8.5	58165	■	SC 3.75 L8.5	58265	■	SCX 3.75 L8.5
57602	■	S 3.75 L10	58166	■	SC 3.75 L10	58266	■	SCX 3.75 L10
57603	■	S 3.75 L11.5	58167	■	SC 3.75 L11.5	58267	■	SCX 3.75 L11.5
57604	■	S 3.75 L13	58168	■	SC 3.75 L13	58268	■	SCX 3.75 L13
			58169	■	SC 3.75 L15	58269	■	SCX 3.75 L15
<b>BEGO Semados® S Implants 4.1</b>			<b>BEGO Semados® SC Implants 4.1</b>			<b>BEGO Semados® SCX Implants 4.1</b>		
57840	■	S 4.1 L7	58170	■	SC 4.1 L7	58270	■	SCX 4.1 L7
57841	■	S 4.1 L8.5	58171	■	SC 4.1 L8.5	58271	■	SCX 4.1 L8.5
57842	■	S 4.1 L10	58172	■	SC 4.1 L10	58272	■	SCX 4.1 L10
57843	■	S 4.1 L11.5	58173	■	SC 4.1 L11.5	58273	■	SCX 4.1 L11.5
57844	■	S 4.1 L13	58174	■	SC 4.1 L13	58274	■	SCX 4.1 L13
			58175	■	SC 4.1 L15	58275	■	SCX 4.1 L15
<b>BEGO Semados® S Implants 4.5</b>			<b>BEGO Semados® SC Implants 4.5</b>			<b>BEGO Semados® SCX Implants 4.5</b>		
57700	■	S 4.5 L7	58176	■	SC 4.5 L7	58276	■	SCX 4.5 L7
57701	■	S 4.5 L8.5	58177	■	SC 4.5 L8.5	58277	■	SCX 4.5 L8.5
57702	■	S 4.5 L10	58178	■	SC 4.5 L10	58278	■	SCX 4.5 L10
57703	■	S 4.5 L11.5	58179	■	SC 4.5 L11.5	58279	■	SCX 4.5 L11.5
57704	■	S 4.5 L13	58180	■	SC 4.5 L13	58280	■	SCX 4.5 L13
			58181	■	SC 4.5 L15	58281	■	SCX 4.5 L15
<b>BEGO Semados® S Implants 5.5</b>			<b>BEGO Semados® SC Implants 5.5</b>			<b>BEGO Semados® SCX Implants 5.5</b>		
57901	■	S 5.5 L8.5	58182	■	SC 5.5 L7	58282	■	SCX 5.5 L7
57902	■	S 5.5 L10	58183	■	SC 5.5 L8.5	58283	■	SCX 5.5 L8.5
57903	■	S 5.5 L11.5	58184	■	SC 5.5 L10	58284	■	SCX 5.5 L10
57904	■	S 5.5 L13	58185	■	SC 5.5 L11.5	58285	■	SCX 5.5 L11.5
			58186	■	SC 5.5 L13	58286	■	SCX 5.5 L13