

**One-stage
overdenture
implant solution**



For various BEGO overdenture solutions

- Machined collar of 3 mm
- Same implant-abutment-interface like BEGO Semados® S/RI implant system
- BEGO CAD/CAM bar restorations
- Sub-Tec Easy-Con abutments
- Sub-Tec Ball attachments

Thread
Self-tapping threads

Surface
High-purity and homogenous
TiPure^{Plus} surface

THE BEGO SEMADOS® S IMPLANTS TG

A high-performance implant system for overdenture implant solutions

BEGO Semados® S implants TG: The last word in modern and cost-effective dental implantology!

- Economic because of short treatment times
- Simple to use in combination with BEGO Semados® S-Line Tray^{Plus}
- Functional
- Quality product manufactured entirely in Germany
- Implant body based on the clinically tested BEGO Semados® S implant
- Prosthetic interface on tissue level

The BEGO Semados® S Implant TG at a glance

- Made from Grade 4 commercially pure titanium
- Diameters: 3.75 · 4.1 mm
- Lengths: 10 · 13 · 15 (with 3 mm machined shoulder)
- TG = Transgingival
- Excellent osseointegration
- High-purity, homogenous TiPure^{Plus} surface
- Rounded apex to protect anatomical structures

The implant-abutment-connection at a glance

- Removable prosthetics (threaded connection)
- Stable connection without micro-movements when subject to functional masticatory loading
- Internal taper for steady seating of the abutment elements for better bacterial resistance
- Customised CAD/CAM bars made from various materials

www.bego-implantology.com

Partners in Progress



My address · Delivery address · Invoice address

I wish to order the following products,
please deliver to the above address.

BEGO Semados® S Implants TG

Art. no.	Colour code	Product	Order quantity
BEGO Semados® S Implants 3.75 TG			
57643	■	S 3.75 TG L 10 (12.5*)	
57644	■	S 3.75 TG L 13 (15.5*)	
57645	■	S 3.75 TG L 15 (17.5*)	
BEGO Semados® S Implants 4.1 TG			
57646	■	S 4.1 TG L 10 (12.5*)	
57647	■	S 4.1 TG L 13 (15.5*)	
57648	■	S 4.1 TG L 15 (17.5*)	

TG = Transgingival

*Overall implant length including 3 mm machined shoulder

Place, date, signature

BEGO Implant Systems GmbH & Co. KG

Wilhelm-Herbst-Str. 1 · 28359 Bremen, Germany

Tel. +49 421 2028-268 · Fax +49 421 2028-265

E-mail info@bego-implantology.com · www.bego-implantology.com